Jet dentaire
Oral Irrigator
Monddouche
Irrigador dental

CONSIGNES RÉGLEMENTAIRES	02
REGULATORY INSTRUCTIONS	08
RICHTLIJNEN	14
CONSIGNAS NORMATIVAS	20

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Please read the user manual after having read the regulatory warnings.

Safety instructions

- Please read the instructions carefully before using the unit and keep them for future reference. If you pass this unit to another person, make sure to include this manual.

DANGER

- To avoid injury resulting from improper use, never use the unit without first having understood how it works and the safety instructions.
- Strictly observe the operating procedure specified in this instruction manual when using the unit.
- To avoid the risk of damage and electric shock, do not disassemble or modify any of the unit's components.
- To avoid burns when cleaning your mouth, do not fill the water tank with water hotter than 40°C.
-  **WARNING! Risk of electric shock! Should the unit fall into water whilst charging, please unplug it first before removing it from the water.**
- If the power cable or plug is damaged, please replace the charger so that it can be used once more.
- Do not heat the unit in a fire. Do not charge, use or dispose of it in a high-temperature environment.
- In the event of a malfunction or obvious damage, please contact our after-sales service.

WARNING

- Only use this unit for its intended purpose, i.e. the one described in this manual.
- To avoid damaging your gums, do not apply excessive force to your gums or teeth with the jet tip.
- Install the unit in a safe place to prevent it from becoming damaged or falling into water.

- To avoid damaging the unit, do not let liquids enter the unit or immerse it in water. Make sure that no foreign objects enter the unit.
- Do not immerse the unit's body or charger in water or any other liquids.
- Do not plug in or unplug the unit with wet hands.
- Make sure that your mains voltage corresponds to that indicated on the rating plate.
- Do not use any unit if the charger, cable or plug is damaged, malfunctions, has been damaged in any way or has fallen into water. Take the unit to an authorised after-sales service centre for inspection, repair or adjustment.
- If the power cable is damaged, it must be replaced by the manufacturer, an authorised service agent or a similarly qualified person, in order to prevent any risks.
- Do not damage, destroy or modify the power cable.
- Do not drag, bend or twist the power cable. Do not place heavy objects on top of it or crush it between other objects.
-  **WARNING! Risk of electric shock! Make sure that the power cable can never get wet or damp when in use. Arrange the power cable in such a way that it cannot be pinched or otherwise damaged.**
- To unplug the unit, do not pull on the power cable; instead grasp the plug and pull it out of the socket.
- Keep the unit and its power cable away from heat-generating objects.
- Unplug the unit before carrying out any maintenance operations.
- Do not install the unit in a high-temperature or humid environment. Avoid exposure to sunlight.
- Do not charge the unit in a high-temperature or humid environment, such as a bathroom or toilet.
- When not in use for a long period of time, please unplug the unit to avoid any fire or other risks.

NOTE:

- Should the unit need to be repaired, please send it to the manufacturer as no parts can be replaced or repaired.
- Please do not use the unit without a jet tip.
- Do not fill the unit with mouthwash; this could cause the unit to corrode.
- This unit may be used by children aged 8 years and over, persons with reduced physical, sensory or mental capacities, and persons with no experience of using or knowledge of the unit, provided they are supervised or receive prior instruction regarding how to safely use the unit and understand the risks involved. Children must not play with the unit. Unit cleaning and maintenance operations must not be performed by unsupervised children.
- Oral irrigators may be used by children aged 8 years and over, persons with reduced physical, sensory or mental capacities, and persons with no experience of using and/or knowledge of the unit, provided they are supervised. Unit cleaning and maintenance operations must not be performed by children unless they are being supervised.
- Children must not play with the unit.
- The unit may only be used with the supplied power adapter.
- The unit's power cable cannot be replaced. Should this cable become damaged, the entire unit must be disposed of.
- If you have had dental surgery in the last two months, or are currently receiving dental treatment, please check with your dentist before using the unit.
- Keep these instructions.
- If the power cable is damaged, it must be replaced by the manufacturer, one of its technicians or a similarly qualified person, in order to prevent any risks.
- Do not let children use the unit as a toy.
-  **WARNING! Risk of electric shock! Never touch the charger if it has fallen into water. Make sure that you unplug the plug from the wall socket first!**

Battery directive

- As the battery is built into the unit, it should only be disposed of by a specialist company.
- Never dispose of this unit with your regular household waste. This product is subject to the provisions of the EU Directive 2012/19/EU.
- Dispose of the unit through an authorised waste collection company or your local recycling centre. Please comply with all the applicable regulations. Please contact your local waste recycling centre if you are unsure.

Disposing of your old device

SEPARATE COLLECTIONS FOR ELECTRONIC AND ELECTRICAL DEVICES

This product must not be disposed of or recycled with other household waste when it reaches the end of its useful life.

The product and its packaging should be collected and recycled separately from other waste.

WHAT SHOULD I DO WITH MY WASTE?

Please dispose of your waste in accordance with the advice below to reduce the impact of your new purchase on the environment.

You are legally responsible for separating the waste associated with this product (and any accessories) so that it can be recycled, including any instructions or/and packaging.

Electric and electronic appliances should not be sorted or recycled together with other household waste.

Keep used products and packaging out of the reach of children, as they can pose a safety hazard. Please retain the instructions and any other notices throughout the service life of the device, and include these documents with the device if it is passed on to third party.

For packaging and other notices, please refer to the recycling instructions on the main packaging, and to the advice available on the Internet.

For the product itself, please refer to the instructions available on the Internet.

If you do not have Internet access, please contact your local authority for advice on how to sort and recycle your waste.

For your children's safety, please store your old devices in a safe place outside your home until they are collected.

SEPARATE COLLECTIONS FOR BATTERIES

The batteries must be changed in accordance with current disposal regulations. Please take the used batteries to a dedicated recycling centre, in order for them to be disposed of in an environmentally friendly manner.

Battery recycling

- This unit contains an internal lithium battery that could explode or release dangerous chemicals.
- To reduce the risk of fire or burns, never disassemble, burn, bury, crush, drill or short-circuit batteries.
- Do not throw the batteries out with household waste.
- To preserve the environment, dispose of batteries in accordance with local regulations.

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